

2-in-1 mobile baby scale and flat scale for toddlers

seca
Precision for health

seca 354

- + Weigh babies lying down and toddlers standing up
- + Compact, robust and lightweight: ideal for mobile use
- + High-quality medical grade baby scale at an economical price
- + Wide range of functions: TARE, HOLD and Breast-Milk-Intake-Function (BMIF)



seca 354

2-in-1 mobile baby scale and flat scale for toddlers



Convertible and robust – baby scale and child scale in one

The weighing tray can be detached from the underframe at the touch of a button and a baby scale becomes a flat scale for toddlers. This allows babies to be weighed lying down and toddlers, to be weighed standing up. The connection is so secure and resilient that you can lift and carry the calibrated seca 354 on the weighing tray.

Capacity	44 lbs / 20 kg
Graduation	0.2 oz < 22 lbs > 0.5 oz / 5 g < 10 kg > 10 g
Dimensions (WxHxD)	21.7 x 6.1 x 13.1 in / 552 x 156 x 332 mm
Tray (WxHxD)	20.7 x 3.1 x 9.8 in / 525 x 80 x 250 mm
Platform (WxHxD)	11 x 2.6 x 9.3 in / 280 x 65 x 235 mm
Net weight	5 lbs / 2.3 kg
Power supply	Batteries
Functions	TARE, Auto-HOLD, Automatic switch-off, lbs/kg switch-over, BMIF
Accessories	+ Measuring mat seca 210 + Measuring board seca 417 + Mobile cart seca 402 + Mobile cart with drawer seca 403 + Backpack seca 409 (for combination with measuring board or measuring mat) + Carrying case seca 413 + Carrying case seca 414 (for combination with measuring board)



Light and mobile – an ideal companion

The detachable weighing tray and the light weight make the compact baby scale an ideal companion, e.g., for midwives on home visits. With a battery capacity of around 10,000 weighings, the scale is ready for use at any location, regardless of the power supply.



Quality that pays off

In addition to the high processing quality and resilience, the baby scale has various practical functions. For example, the breast milk intake function (BMIF), which determines the amount of milk absorbed by the baby with the baby scale calculating the difference between the weight before and after food intake. This makes the scale ideal for parents who want to control the development of their child after birth at home.

We reserve the right to make modifications.

seca north america | west
13601 Benson Avenue
Chino, CA 91710
USA

seca north america | east
1439 Stuart Engals Blvd, Suite 200
Mount Pleasant, SC 29464
USA

T +1 800 542 7322
F +1 888 705 7397
E info.us@seca.com

Precision for health
seca.com